

LUS Image Classification of Covid-19 and Pneumonia

Abhijith C

242CS003

Under the guidance of **Dr. Jeny Rajan**



March 24, 2025

Introduction

Lung ultrasound (LUS) images provide key diagnostic information by capturing different artifacts and structures in the lungs. Key elements in LUS images include:

- **Pleural Line** – The bright, horizontal line at the top of the image, representing the lung's surface.
- **A-Lines** – Repetitive horizontal artifacts indicating normal aerated lungs.
- **B-Lines** – Vertical, hyperechoic (bright) lines extending to the bottom, associated with lung pathologies like pneumonia and pulmonary edema.

LUS Images



LUS Differences: COVID-19 vs. Pneumonia

Feature	COVID-19	Pneumonia
B-Lines	Scattered, widespread	Localized in one area
Pleural Line	Thick and uneven	Mostly normal
Lung Involvement	Both lungs (bilateral)	One lung (unilateral)

Dataset

Dataset Overview

Label	Count
COVID-19	524
Pneumonia	463

Table 1: Dataset Label Distribution

Data Loading Process:

- Resize images to uniform dimensions (224×224).
- Convert labels into numerical format (COVID = 0, Pneumonia = 1).

Data Splitting:

- Train (80%) + Validation (10%) + Test (10%) split.

Data Augmentation Steps

- **Geometric Transformations:**
 - Rotation: Randomly rotates images by up to $\pm 20^\circ$.
 - Width Shift: Shifts images horizontally by $\pm 10\%$.
 - Height Shift: Shifts images vertically by $\pm 10\%$.
 - Shear: Applies shearing transformation up to 0.2 radians.
 - Zoom: Random zooming in/out up to $\pm 20\%$.
 - Horizontal Flip: Randomly flips images left-right.
- **Fill Mode:** Uses nearest pixel values to fill missing areas.

Experiments Performed

Model Architecture:

- **Base Model:** Pretrained ResNet50 model with weights trained on ImageNet
- **Custom Layers:**
 - Flatten Layer
 - Fully Connected Layer (512 units, ReLU activation)
 - Output Layer (1 unit, Sigmoid activation)
- **Loss Function:** Binary Cross-Entropy
- **Optimizer:** Adam (Learning Rate = 0.0001)

ResNet50




Figure: Accuracy Curve

ResNet50




Figure: Loss Curve

ResNet50




Figure: Confusion Matrix

ResNet50

Final Test Accuracy: 1.0000

Final Test Loss: 0.0003

Classification Report for ResNet50:

	precision	recall	f1-score	support
covid	1.00	1.00	1.00	50
pneumonia	1.00	1.00	1.00	49
accuracy			1.00	99
macro avg	1.00	1.00	1.00	99
weighted avg	1.00	1.00	1.00	99

Model Architecture:

- **Base Model:** Pretrained VGG16 model with weights trained on ImageNet
- **Custom Layers:**
 - Flatten Layer
 - Fully Connected Layer (512 units, ReLU activation)
 - Output Layer (1 unit, Sigmoid activation)
- **Loss Function:** Binary Cross-Entropy
- **Optimizer:** Adam (Learning Rate = 0.0001)

VGG16




Figure: Accuracy Curve




Figure: Loss Curve




Figure: Confusion Matrix

Final Test Accuracy: 1.0000

Final Test Loss: 0.0033

Classification Report for VGG16:

	precision	recall	f1-score	support
covid	1.00	1.00	1.00	50
pneumonia	1.00	1.00	1.00	49
accuracy			1.00	99
macro avg	1.00	1.00	1.00	99
weighted avg	1.00	1.00	1.00	99

K-Fold Cross-Validation for ResNet50

Key Steps:

- **Splitting Data:** The dataset is divided into 5 folds.
- **Training & Validation:** Each fold is used once as validation while the remaining 4 folds are used for training.
- **Performance Evaluation:** Validation accuracy and loss are computed for each fold.
- **Final Model Training:** After 5 iterations, a final model is trained on the full dataset.
- **Test Set Evaluation:** The final model is evaluated on the test set.

K-Fold Cross-Validation for ResNet50

Fold	Validation Accuracy	Validation Loss
Fold 1	0.9430	0.6063
Fold 2	0.9937	0.0074
Fold 3	0.9937	0.0267
Fold 4	0.9494	1.8617
Fold 5	0.8280	40.4298

Average Validation Accuracy across 5 folds: 0.9416
Average Validation Loss across 5 folds: 8.5864
Average Training Time per fold: 117.22 seconds

K-Fold Cross-Validation for ResNet50




Figure: Confusion Matrix

K-Fold Cross-Validation for ResNet50

Final Test Accuracy: 0.9495

Final Test Loss: 6.3086

Classification Report for KFold_Model:

	precision	recall	f1-score	support
covid	1.00	0.90	0.95	105
pneumonia	0.90	1.00	0.95	93
accuracy			0.95	198
macro avg	0.95	0.95	0.95	198
weighted avg	0.95	0.95	0.95	198

- **Original Class Distribution:**

- Class 0 (Covid): 524 samples
- Class 1 (Pneumonia): 463 samples

- **Imbalanced Class Distribution:**

- Class 0 (Covid): 524 samples
- Class 1 (Pneumonia): 46 samples

ResNet50 with Imbalanced Dataset




Figure: Accuracy Curve

ResNet50 with Imbalanced Dataset



Figure: Loss Curve

ResNet50 with Imbalanced Dataset




Figure: Confusion Matrix

ResNet50 with Imbalanced Dataset

Final Test Accuracy: 1.0000

Final Test Loss: 0.0057

Classification Report for ResNet50:

	precision	recall	f1-score	support
covid	1.00	1.00	1.00	50
pneumonia	1.00	1.00	1.00	7
accuracy			1.00	57
macro avg	1.00	1.00	1.00	57
weighted avg	1.00	1.00	1.00	57

Introduction to Vision Transformers (ViTs)

- Originally designed for **Natural Language Processing (NLP)**, used in models like **GPT** and **BERT**.
- **First applied to computer vision** in 2020 for **image classification and segmentation** by Dosovitskiy et al. [3].
- The **core idea**:
 - Instead of convolutional filters, **images are divided into patches**.
 - Each patch is treated as a **token**, similar to words in NLP.
 - A standard **Transformer encoder** is used to process these tokens.
- Great potential in Ultrasound Image Analysis [2].

ViT Architecture




Figure: Model Overview of Vision Transformers

Introduction to Swin Transformer

- Swin Transformer was proposed by Microsoft Research in 2021 [1].
- The **core idea**:
 - Uses a **shifted window attention mechanism** instead of global self-attention.
 - Reduces computational complexity compared to ViT.
- Captures Fine Details with Local & Global Features.
- Scales well and processes high resolution images faster.
- Used in COVID-19 Pneumonia Assessment in Lung Ultrasound Images [3].

Swin Tiny Transformer

- **Base Model:** Swin Tiny Transformer pretrained on ImageNet-1K
- **Input Size:** $224 \times 224 \times 3$
- **Frozen Layers:** All except classification head
- **New Classification Head:**
 - Linear (num_features \rightarrow 256)
 - ReLU Activation
 - Dropout (0.5)
 - Linear (256 \rightarrow 2 classes)
- **Loss Function:** CrossEntropyLoss
- **Optimizer:** Adam (lr= 1×10^{-4} , weight decay= 1×10^{-5})
- **Scheduler:** StepLR (Step size = 5, Gamma = 0.5)

Swin Tiny Transformer




Figure: Accuracy Curve

Swin Tiny Transformer





Figure: Loss Curve

Swin Tiny Transformer



Test Accuracy: 98.48%




Test Loss: 0.0783

Conclusion




Performance Metrics for Different Models

Model	Train Accuracy	Test Accuracy	Train Loss	Test Loss
ResNet50	99.69	100	0.0413	0.003
VGG16	97.59	100	0.1192	0.033
ResNet50 K-Fold	99.60	94.95	0.0123	6.3086
ResNet50 Imbalanced	97.79	100	0.2008	0.0057
Swin Transformer	96.60	98.48	0.1477	0.0783

References

-  Xing, W., et al. "Automatic detection of A-line in lung ultrasound using deep learning." Medical Physics 50.1 (2023): 330-343.
-  Baloesu, C., et al. "Automated lung ultrasound B-line assessment using deep learning." IEEE T-UFFC 67.11 (2020): 2312-2320.
-  Dosovitskiy, A., et al. "An image is worth 16x16 words: Transformers for image recognition at scale." arXiv:2010.11929 (2020).

References

-  Liu, Z., et al. "Swin Transformer: Hierarchical ViT using shifted windows." ICCV (2021).
-  Vafaezadeh, M., Behnam, H., & Gifani, P. "Ultrasound image analysis with vision transformers." Diagnostics 14.5 (2024): 542.
-  Fiorentino, M. C., et al. "ViT approaches for COVID-19 pneumonia assessment in lung ultrasound." IEEE MetroXRAINE (2024).

Thank You!